

PATELLA FIXATION BANDAGE

PATELLA CORRECTION BANDAGE

KNEE BANDAGES

StabiloGen® ...the classical bandage!



BORT StabiloGen®

No. 114 420



BORT Knee Bandage with Patella Recess

No. 114 150

- Patella fixation bandage
- Integrated silicone ring
- 2 splints each for improved lateral stability
- Strong dual tension quality
- With padding for the hollow of the knee to provide optimal comfort
- Micro fibre material

Colors: beige, silvergrey, blue

Indication: Patella instability (luxation tendency), articular effusion, swellings, arthrosis and arthritis, post-operative irritations, ligament instability

- Patella ring of silicone
- 1 splint each side
- Strong dual tension quality
- With pressure-reducing edges
- Micro fibre material

Colors: beige, silvergrey

Indication: Chondropathia patellae, instabilities of knee ligaments of different nature, meniscopathia, post-operative, bursitis

Sizes	Circumference of knee
small	up to 34 cm
medium	34 – 37 cm
large	37 – 40 cm
x-large	40 – 44 cm
xx-large	44 – 48 cm

Asymmetric ...the correction bandage!



BORT Asymmetric

No. 114 700

- Individual guiding of patella towards medial and distal by elastic reins
- Anatomically form-knitted
- Micro fibre material
- Patella silicon ring
- 1 splint each side
- Pressure-reducing zones in the peripheral areas and the hollow of the knee

Color: silvergrey

Indication: Patella laterality, luxation tendency, condition after lateral release operation. Patella instability, meniscopathies, chondropathia patellae.

Sizes	Circumference of knee
small	up to 34 cm
medium	34 – 37 cm
large	37 – 40 cm
x-large	40 – 44 cm
Please indicate left or right	

STABILOGen® Zip



NEW
THE KNEE BANDAGE WITH
ZIP


ZAG-Knee/GB-01/05 mm

BORT GmbH
P.O.Box 1330
D-71367 Weinstadt-Benzach
www.bort.de
Phone **49 (0) 71 51-9 92 00-0



StabiloGen® Zip

... the comfortable with zip!



NOVELTY
with zip fastener
for easy putting on
and taking off.



Comfortable putting on:
Attach bandage to calf and close zip fastener.
Pull-up bandage in closed condition.

Easy taking off:
Now very simple without low bending down. Open zip fastener completely – finished.

BORT StabiloGen® Zip

No. 114 390

- Patella fixation bandage
- With zip fastener for easy putting-on and taking off
- Integrated silicone ring
- 1 splint each for improved lateral stability
- Strong dual tension quality
- With padding for the hollow of the knee to provide optimal comfort
- Micro fibre material

Colors: beige, silvergry, blue

Indication:
Patella instability (luxation tendency), articular effusions, swellings, arthrosis and arthritis, post-operative irritations, ligament instability

Sizes	Circumference of knee
small	up to 34 cm
medium	34 – 37 cm
large	37 – 40 cm
x-large	40 – 44 cm
xx-large	44 – 48 cm



Perfect fit
The StabiloGen® Zip
also during sports activities,
e.g. during semi-marathon.